

Appendix B2

**Draft ICCVAM Recommended Protocol for Future Studies Using the Mono Mac 6
(MM6)/Interleukin-6 (IL-6) Test Method**

[This Page Intentionally Left Blank]

1 **Draft ICCVAM Recommended Protocol for Future Studies Using the Mono Mac 6**
2 **(MM6)/Interleukin-6 (IL-6) Test Method**

3

4

5
6 **PREFACE**

7 This proposed protocol for the detection of pyrogenicity is based on information obtained
8 from 1) The European Centre for the Validation of Alternative Methods (ECVAM) MM6/IL-
9 6 Background Review Document (BRD) presented in Appendix A of the draft Interagency
10 Coordinating Committee on the Validation of Alternative Methods (ICCVAM) Pyrogenicity
11 Test Method BRD, which includes ECVAM Standard Operating Procedures (SOPs) for the
12 MM6/IL-6 test method, and 2) Information provided to the National Toxicology Program
13 (NTP) Interagency Center for the Evaluation of Alternative Toxicological Methods
14 (NICEATM) by Dr. Thomas Hartung, Head of ECVAM. The ECVAM SOPs are based on
15 the method published by Taktak et al. (1991). A table of comparison between the draft
16 ICCVAM recommended protocol and the ECVAM SOPs is provided in **Table 1**. Future
17 studies using the MM6/IL-6 test method may include further characterization of the
18 usefulness or limitations of the assay for regulatory decision-making. Users should be aware
19 that the proposed test method protocol might be revised based on additional optimization
20 and/or validation studies. ICCVAM recommends that test method users routinely consult the
21 ICCVAM/NICEATM website (<http://iccvam.niehs.nih.gov/>) to ensure that the most current
22 test method protocol is used.

22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38

[This Page Intentionally Left Blank]

39 **Table 1 Comparison of Draft ICCVAM Recommended Test Method Protocol with the ECVAM SOP for the Mono Mac**
 40 **6 (MM6)/Interleukin-6 (IL-6) Pyrogen Test Method**

Protocol Component	Draft ICCVAM Protocol	ECVAM SOP	ECVAM Validation SOP
Test Substance	Test neat or at minimal dilution that produces no interference	Same as ICCVAM protocol	Test at MVD
Incubation Plate (number of control or test groups at n=4 each)	NSC (1)	Same as ICCVAM protocol	Same as ICCVAM protocol
	EC (5)	Same as ICCVAM protocol	Same as ICCVAM protocol
	TS (14)	Same as ICCVAM protocol	TS (2) x EC (5) spikes
	PPC ¹ (0)	Same as ICCVAM protocol	PPC (2)
	NPC ¹ (0)	Same as ICCVAM protocol	NPC (2)
	PC (0)	Same as ICCVAM protocol	PC (1)
	NC (0)	Same as ICCVAM protocol	NC (1)
ELISA Plate	Includes seven point IL-6 SC and blank in duplicate	Same as ICCVAM protocol	Same as ICCVAM protocol
Decision Criteria for Interference	0.5 x Median OD ₄₅₀ ² of 1 EU/mL EC <2x Median OD ₄₅₀ of 1 EU/mL EC	Same as ICCVAM protocol	Not applicable (tested at MVD)
Assay Acceptability Criteria	Mean OD ₄₅₀ ² of PPC is 50% to 200% of 0.5 EU/mL EC	Same as ICCVAM protocol	Same as ICCVAM protocol
	Mean OD ₄₅₀ of NSC ≤0.15	Same as ICCVAM protocol	Mean OD ₄₅₀ of NSC ≤0.20
	Not included	Outliers rejected using Dixon's text (p≥0.5)	Outliers rejected using Dixon's text (p≥0.5)
	Not included	Mean OD ₄₅₀ EC > NSC (2SD with n-1 weighting)	Mean OD ₄₅₀ PC > LOQ ³
	Quadratic function of IL-6 SC r ² ≥0.95	Same as ICCVAM protocol	Same as ICCVAM protocol
	EC SC produces OD ₄₅₀ values that ascend in a sigmoidal concentration response	Same as ICCVAM protocol	Mean OD ₄₅₀ of each EC > Mean OD ₄₅₀ of next lower EC concentration (minimum of 4 data points needed for valid SC)
Decision Criteria for Pyrogenicity	Mean OD ₄₅₀ ² of TS > Mean OD ₄₅₀ of 0.5 EU/mL EC	EC SC data transformed to 4-parameter logistical model by an in-house program and the linear mean square is calculated. TS pyrogen content is compared with the ELC ⁴ using confidence limits for significance. Dixon's test is used to reject outliers.	EC SC data transformed to 4-parameter logistical model by an in-house program and the linear mean square is calculated. TS pyrogen content is compared with the ELC ⁴ using confidence limits for significance. Dixon's test is used to reject outliers.

41 Abbreviations: EC = Endotoxin control; ELC = Endotoxin Limit Concentration; LOQ = Limit of Quantification; MVD = Maximum Valid Dilution; NC =
42 Negative Control; NPC = Negative Product Control; NSC = Normal saline control; PC = Positive control; PPC = Positive Product Control; SC = Standard curve;
43 SD = Standard Deviation; TS = Test substance
44 ¹ PPC and NPC are evaluated during the interference test.
45 ² Median or mean OD₄₅₀ values are corrected (i.e., reference filter reading, if applicable, and NSC are subtracted).
46 ³ LOQ is the mean OD₄₅₀ of the NSC + 10xSD mean OD₄₅₀ of the NSC.
47 ⁴ Where unknown, the ELC is calculated (e.g., Based on a rabbit sensitivity of 5 EU/kg, for a product injected at 10 mL/kg, the detection limit is 5 EU/10 mL/kg
48 or 0.5 EU/mL/kg, or an ELC of 0.5 EU/mL).
49
50

51 **1.0 PURPOSE AND APPLICABILITY**

52 The purpose of this protocol is to describe the procedures used to evaluate the presence of a
53 pyrogen (i.e., Gram-negative endotoxin) in a test substance. The presence of Gram-negative
54 endotoxin is detected by its ability to induce cytokine IL-6 release from MM6 cells. The
55 quantity of IL-6 released is obtained using an enzyme-linked immunosorbent assay (ELISA)
56 that includes monoclonal or polyclonal antibodies specific for IL-6. Release of this cytokine
57 is measured by incubation of the MM6 cells with test substances or controls (i.e., positive,
58 negative). The amount of pyrogen present is determined by comparing the values of
59 endotoxin equivalents produced by MM6 cells exposed to the test substance to those exposed
60 to an internationally-harmonized Reference Standard Endotoxin (RSE)¹ or an equivalent
61 standard expressed in Endotoxin Units (EU)/mL. Based on a rabbit threshold pyrogen dose of
62 0.5 EU/mL, which was established in a retrospective evaluation of rabbit pyrogen test (RPT)
63 data, a test substance is considered pyrogenic if it induces a level of IL-6 release equal to or
64 greater than 0.5 EU/mL.

65 The focus of this protocol is on the use of the MM6/IL-6 test method specifically for the
66 detection of Gram-negative endotoxin in parenteral pharmaceuticals. The relevance and
67 reliability for non-endotoxin pyrogens (e.g., lipoteichoic acid) has not been demonstrated in a
68 formal validation study.

69 **2.0 SAFETY AND OPERATING PROCEDURES**

70 All procedures for procurement of eligible blood donors and blood donations should follow
71 the regulations and procedures set forth by institutional guidelines for utilization of human
72 substances, which include but are not limited to blood, tissues, and tissue fluids. Standard
73 laboratory precautions are recommended, including the use of laboratory coats, eye

¹ RSEs are internationally-harmonized reference standards (e.g., WHO *E. coli* Lipopolysaccharide [LPS] 94/580 [0113:H10:K-]; USP RSE Lot G3E069; FDA Lot EC-6). Equivalent endotoxins include commercially available *E. coli*-derived LPS Control Standard Endotoxin (CSE) or other *E. coli* LPS preparations that have been calibrated with an appropriate RSE.

74 protection, and gloves. If necessary, additional precautions required for specific study
75 substances or hazardous chemicals will be identified in the Material Safety Data Sheet
76 (MSDS).

77 The stop solution used in the ELISA kit is acidic and corrosive and should be handled with
78 the proper personal protective devices. If this reagent comes into contact with skin or eyes,
79 wash thoroughly with water. Seek medical attention, if necessary.

80 Tetramethylbenzidine (TMB) solution contains a hydrogen peroxide substrate and 3, 3', 5,
81 5'-TMB. This reagent is a strong oxidizing agent and a suspected mutagen. Appropriate
82 personal protection should be used to prevent bodily contact.

83 Bacterial endotoxin is a toxic agent (i.e., can induce sepsis, shock, vascular damage,
84 antigenic response) and should be handled with care. Skin cuts should be covered and
85 appropriate personal protective devices should be worn. In case of contact with endotoxin,
86 immediately flush eyes or skin with water for at least 15 min. If inhaled, remove the affected
87 individual from the area and provide oxygen and/or artificial respiration as needed. Skin
88 absorption, ingestion, or inhalation may produce fever, headache, and hypotension.

89 **3.0 MATERIALS, EQUIPMENT AND SUPPLIES**

90 **3.1 Source of Cells**

91 The MM6 cell line is a human monocytic cell line originally described by Professor H.W.L.
92 Ziegler-Heitbrock (Institute for Immunology, University of Munich, Munich, Germany). A
93 Master Cell Bank and a Working Cell Bank have been established at the National Institute
94 for Biological Standards and Control (NIBSC), from which the MM6 cells can be purchased.

95 **3.2 Equipment and Supplies**

96 For all steps in the protocol, excluding the ELISA procedure, the materials that will be in
97 close contact with samples (e.g., pipet tips, containers, solutions) should be sterile and free
98 from detectable pyrogens.

- 99 3.2.1 Utilization of MM6 cells
- 100 3.2.1.1 Equipment
- 101 • Centrifuge
- 102 • Hemacytometer
- 103 • Hood; Bio-safety, laminar flow (recommended)
- 104 • Incubator; cell culture ($37\pm 1^{\circ}\text{C} + 5\% \text{CO}_2$)
- 105 • Inverted Microscope
- 106 • pH meter
- 107 • Pipetter; multichannel (8- or 12-channel)
- 108 • Pipettors; single-channel adjustable (20, 200, and 1000 μL)
- 109 • Repeating pipetter
- 110 • Vortex mixer
- 111 • Water bath
- 112 3.2.1.2 *Consumables*
- 113 • Centrifuge tubes; nonpyrogenic, polystyrene (15 and 50 mL)
- 114 • Cryotubes; screw-cap (2 mL)
- 115 • Combitips; repeating pipetter (1.0 and 2.5 mL)
- 116 • Filters; sterile, 0.22 μm
- 117 • Flasks; tissue culture
- 118 • Phosphate buffered saline (PBS); sterile
- 119 • Pipets; sterile, glass
- 120 • Plates; microtiter, nonpyrogenic, 96-well, polystyrene, tissue culture
- 121 • Reaction tubes; polystyrene (1.5 mL)
- 122 • Reservoirs; fluid

- 123 • RPMI-1640 cell culture medium; supplemented as described in **Section 4.3** to
124 yield RPMI-Complete (RPMI-C)
- 125 • Tips; pipetter, sterile, pyrogen-free (20 and 200 μ L)
- 126 • Tubes; polystyrene

127 3.2.2 ELISA

128 3.2.2.1 *Equipment*

- 129 • Microplate mixer
- 130 • Microplate reader (450 nm with an optional reference filter in the range of
131 540-590 nm)
- 132 • Microplate washer (optional)
- 133 • Multichannel pipetter

134 3.2.2.2 *Consumables*

- 135 • Container; storage, plastic
- 136 • Deionized water; nonsterile
- 137 • Plates; microtiter, nonpyrogenic, 96-well, polystyrene, tissue culture
- 138 • Pyrogen-free water (PFW)
- 139 • Reservoirs; fluid
- 140 • Tips; pipetter, nonsterile
- 141 • Tubes; polystyrene (12 mL)

142 3.2.2.3 *ELISA Kit*

143 An ELISA that measures IL-6 release from MM6 cells is used. A variety of IL-6 ELISA kits
144 are commercially available and the IL-6 ELISA procedure outlined in this protocol is
145 intended to serve as an example for using an ELISA kit. If the user prefers to prepare an in-
146 house ELISA, then additional reagents would be required. The IL-6 ELISA should be
147 calibrated using an IL-6 international reference standard (e.g., WHO 89/548) prior to use.
148 The IL-6 cytokine assay kits do not provide the RSE or endotoxin equivalent; therefore, this

149 reagent must be purchased separately. Results obtained using these products are subject to
150 the assay acceptability and decision criteria described in **Sections 8.0** and **9.0**. IL-6 ELISA
151 kit components may include the following:

- 152 • ELISA plates coated with anti-human IL-6 capture antibody; monoclonal or
153 polyclonal
- 154 • Buffered wash solution
- 155 • Dilution buffer
- 156 • Enzyme-labeled detection antibody
- 157 • Human IL-6 reference standard
- 158 • Pyrogen-free saline (PFS)
- 159 • Stop solution
- 160 • TMB/substrate solution

161 **3.3 Chemicals**

- 162 • Endotoxin (e.g., WHO *E. coli* LPS 2nd International Standard 94/580; USP
163 RSE *E. coli* LPS Lot G3E069; USP RSE *E. coli* Lot G; FDA *E. coli* Lot EC6)

164 **3.4 Solutions**

- 165 • RPMI-C cell culture medium

166 **4.0 ASSAY PREPARATION**

167 All test substances, endotoxin, and endotoxin-spiked solutions should be stored at 4°C.

168 **4.1 Endotoxin Standard Curve**

169 An internationally harmonized RSE or equivalent is used to generate the endotoxin standard
170 curve. The use of any other *E. coli* LPS requires calibration against a RSE using the
171 MM6/IL-6 test method.

172 A standard endotoxin curve consisting of a Normal Saline Control (NSC) and five RSE
173 concentrations (0.125, 0.25, 0.50, 1.0, and 2.0 EU/mL) are included in the incubation step

174 (refer to **Table 4-1**) and then transferred to the ELISA plate. To prepare the endotoxin
 175 standard curve, first obtain a 2000 EU/mL stock solution by addition of PFW to the
 176 lyophilized content of the stock vial by following the instructions provided by the
 177 manufacturer (e.g., for a vial containing 10,000 EU, 5 mL of PFW is added). To reconstitute
 178 the endotoxin, the stock vial should be vortexed vigorously for at least 30 min or sonicated in
 179 a bath sonicator for 5 min. The stock solution is stable for 14 days when stored at 2 to 8°C.
 180 An endotoxin standard curve is prepared by making serial dilutions of the stock solution in
 181 PFS as described in **Table 4-1**.

182 **Table 4-1 Preparation of Endotoxin Standard Curve**

Stock Endotoxin EU/mL	µL of Stock Endotoxin	µL of PFS	Endotoxin Concentration EU/mL
2000 ^{1,2}	20	1980	20 ³
20	100	900	2.0
2.0	500	500	1.0
1.0	500	500	0.50
0.50	500	500	0.25
0.25	500	500	0.125
0	0	1000	0

183 Abbreviations: EU = Endotoxin units; PFS = Pyrogen-free saline

184 Each stock tube should be vortexed vigorously prior to its use to make the subsequent dilution.

185 ¹ A 2000 EU/mL stock solution of endotoxin is prepared according to the manufacturer's instructions.

186 ² The stock solution of USP RSE may be stored in aliquots and kept at -20°C for up to 6 months. Do not store
 187 the endotoxin at -80°C.

188 ³ This concentration is not used in the assay.

189

190 **4.2 Test Substances**

191 Liquid test substances should be tested neat or, if interference is detected (see **Section 4.2.1**),
 192 diluted in PFS. Solid test substances should be prepared as solutions in PFS or, if insoluble in
 193 saline, dissolved in DMSO then diluted up to 0.5% (v/v) with PFS provided that this
 194 concentration does not interfere with the assay. The test substances should be vortexed
 195 vigorously for at least 30 min or sonicated in a bath sonicator for 5 min.

196 **4.2.1 Interference Testing**

197 Interference testing must be carried out on any test sample for which no interference
 198 information is available. The purpose of the interference test is to determine the lowest
 199 dilution (i.e., highest concentration) of a test substance from which an endotoxin spike can be

200 detected (i.e., based on the decision criteria described in **Section 4.2.1.2**). However, to ensure
201 a valid test, a test substance should not be diluted beyond its Maximum Valid Dilution
202 (MVD).

203 For many marketed products, values for the MVD and the Endotoxin Limit Concentration
204 (ELC) are published in the U.S. Pharmacopeia, the European Pharmacopoeia, and/or Food
205 and Drug Administration (FDA) guidelines. However, if one or both of these values are not
206 available, then calculation of the MVD is dependent on the ELC (see **Section 12.3**). If
207 unknown, the ELC can be approximated by dividing the maximum hourly dose of the
208 product by the hourly dose received per patient. For example, if a product is used at an
209 hourly dose of 100 mg per patient, then the ELC would be 350 EU/100 mg or 3.5 EU/mg.

210 4.2.1.1 *Reference Endotoxin for Spiking Test Substances*

211 The WHO-LPS 94/580 [*E. coli* O113:H10:K-] or equivalent international RSE (e.g., USP
212 G3E069, FDA EC-6) is recommended for preparation of the endotoxin control (EC). If a
213 different *E. coli* LPS is used and the potency relative to the RSE is not provided, then each
214 lot must be calibrated against the RSE in the MM6/IL-6 test method. For interference testing,
215 an endotoxin standard curve (see **Section 4.1**) should be included on each plate.

216 4.2.1.2 *Spiking Test Substances with Endotoxin*

217 For interference testing, non-spiked and endotoxin-spiked test substances are prepared in
218 microplate wells (n=4 replicates) and an *in vitro* pyrogen test is performed. Either RPMI-C or
219 a fixed concentration (a concentration selected from the middle of the EC standard curve) of
220 the RSE (i.e., 1 EU/mL) in RPMI-C is added to the test substance in serial two-fold dilutions.
221 An illustrative example of endotoxin spiking solutions is shown in **Table 4-2**. For non-spiked
222 solutions, 150 μ L of RPMI-C and 50 μ L of the test substance (neat or at serial dilution) are
223 added to a well. Then, 50 μ L of MM6 cells are added and the well contents mixed are mixed.
224 Endotoxin-spiked solutions are prepared by adding 100 μ L of RPMI-C, 50 μ L of the test
225 substance (neat or at serial dilution), and 50 μ L of an endotoxin spike solution (1.0 EU/mL).
226 MM6 cells (50 μ L) are then added and the well contents are mixed (see example presented in
227 **Table 4-2**).

228

229 **Table 4-2 Preparation of Endotoxin Spiked and Non-Spiked Solutions for**
 230 **Determination of Test Substance Interference in the Incubation and**
 231 **ELISA Test Systems**

Sample Addition	Spiked	Non-spiked
	μL/well ¹	
RPMI-C (containing 2% FCS)	100	150
Endotoxin spike solution ²	50	0
Test substance (neat and each serial dilution)	50	50
MM6 cells ³	50	50
Total ⁴	250	250

232 Abbreviations: MM6 cells = Mono Mac 6 cells

233 ¹ n=4 replicates each

234 ² Endotoxin concentration is 1.0 EU/mL in RPMI-C.

235 ³ MM6 cells are resuspended in RPMI-C (2.5 x 10⁶ cells/mL).

236 ⁴ A total volume of 250 μL per well is used for the incubation.

237

238 The lowest dilution of the test substance that yields an endotoxin spike recovery of 50% to
 239 200% in the pyrogen test is determined. The optical density (OD) values of the endotoxin-
 240 spiked and non-spiked test substances are calibrated against the endotoxin calibration curve.
 241 The resulting EU value of the non-spiked test substance is subtracted from the corresponding
 242 EU value of the endotoxin-spiked test substance at each dilution. The % recovery for each
 243 sample dilution is then determined from the endotoxin spike solution concentration set to
 244 100%. For example, consider the following interference test results in **Table 4-3**:

245 **Table 4-3 Example of Interference Data Used to Determine Sample Dilution for**
 246 **Assay**

Sample Dilution	% Recovery of Endotoxin Control
None	25
1:2	49
1:4	90
1:8	110

247

248 Based on these results, the dilution of the test substance used in the *in vitro* pyrogen test
 249 would be 1:4 (i.e., the lowest dilution between 50% and 200% of the 1.0 EU/mL EC).

250 4.2.2 Interference with ELISA

251 If the data obtained from the experiment in **Section 4.2.1** suggests the presence of
 252 interference, then a subsequent experiment similar to that described in **Section 4.2.1** would

253 need to be performed to confirm that the test substance(s) does not directly interfere with the
254 ELISA. For this experiment, an ELISA would be performed in the absence of MM6 cells.

255 4.3 Cell Culture Medium

256 MM6 cells are maintained in RPMI-C. For use in the ELISA procedure, the concentration of
257 fetal calf serum (FCS) is reduced to 2%. Each medium is prepared and stored as described by
258 the manufacturer.

259 4.3.1 RPMI-C

- 260 • Bovine insulin; 0.23 IU/mL
- 261 • FCS; heat-inactivated (50 mL or a 10% final concentration)
- 262 • HEPES buffer; 20 mM
- 263 • L-Glutamine; 2 mM
- 264 • MEM non-essential amino acids; 0.1 mM
- 265 • Oxaloacetic acid; 1 mM
- 266 • Penicillin/streptomycin (10,000 IU/mL penicillin, 10 mg/mL streptomycin)
- 267 • RPMI-1640 medium (500 mL)
- 268 • Sodium pyruvate; 1 mM

269 4.3.2 Starting a Culture of MM6 Cells

270 To initiate a culture of MM6 cells, remove a vial of the primary stock from liquid nitrogen.
271 Thaw the vial on ice. Gently mix and transfer the cells to a 50 mL centrifuge tube and add 10
272 mL of RPMI-C. Centrifuge at 100 x g for 5 min at room temperature (RT). Remove the
273 supernatant and resuspend the cells in ice-cold RPMI-C. Centrifuge at 100 x g for 5 min at
274 RT. Remove the supernatant and resuspend the MM6 cells in 2 mL of RPMI-C. Add 8 mL of
275 RPMI-C to a 25 cm² tissue culture flask and transfer the cell suspension to the flask. Place
276 the flasks in a cell culture incubator and maintain the cells at 37±1°C + 5% CO₂.

277 4.3.3 Propagation of MM6 Cells

278 To propagate the MM6 cells, transfer the MM6 cells to a 50 mL tube. Centrifuge at 100 x g
279 for 8 min at RT. Remove the supernatant, resuspend the cell pellet in 4 mL of RPMI-C, and
280 gently pipet up and down to mix. Transfer an aliquot of the cell suspension to new tissue
281 culture flasks and add fresh RPMI-C to obtain a final concentration of 2×10^5 cells/mL.
282 Place the flasks in a cell culture incubator and maintain the cells at $37 \pm 1^\circ\text{C} + 5\% \text{CO}_2$.

283 4.3.4 Preparation of a MM6 Cell Bank

284 To initiate a bank of MM6 cells, centrifuge the cell culture(s) at 100 x g for 8 min at 4°C .
285 Remove the supernatant and resuspend the cells in FBS at 4°C . Adjust the cell concentration
286 to 4×10^6 cells/mL and store on ice for 10 min. Add an equal volume of ice-cold FBS
287 containing 10% dimethylsulfoxide (DMSO) drop-wise to the cell suspension (final
288 concentration is 2×10^6 cells/mL with 5% DMSO). Transfer the cell suspension to sterile,
289 pyrogen-free cryotubes (1 mL/tube). Place the tubes in a well-insulated polystyrene box and
290 store at -80°C or below for greater than 48 hr and then transfer to a liquid nitrogen container.

291 **5.0 CONTROLS**

292 **5.1 Negative Control**

293 A negative control (e.g., RPMI-C) is included in each experiment in order to detect
294 nonspecific changes in the test system, as well as to provide a baseline for the assay
295 endpoints.

296 **5.2 Solvent Control**

297 Solvent controls are recommended to demonstrate that the solvent is not interfering with the
298 test system when solvents other than PFS are used to dissolve test substances.

299 **5.3 Positive Control**

300 An internationally standardized EC (e.g., WHO 94/580; USP G3E069; 0.5 EU/mL) is
301 included in each experiment to verify that an appropriate response is induced.

302 **5.4 Benchmark Control**

303 Benchmark controls may be used to demonstrate that the test method is functioning properly,
304 or to evaluate the relative pyrogenic potential of chemicals (e.g., parenteral pharmaceuticals,
305 medical device eluates) of a specific class or a specific range of responses, or for evaluating
306 the relative pyrogenic potential of a test substance. Appropriate benchmark controls should
307 have the following properties:

- 308 • consistent and reliable source(s) for the chemicals (e.g., parenteral
309 pharmaceuticals, medical device eluates)
- 310 • structural and functional similarities to the class of the substance being tested
- 311 • known physical/chemical characteristics
- 312 • supporting data on known effects in animal models
- 313 • known potency in the range of response

314 **5.5 Positive Product Control (PPC)**

315 The PPC is a test substance diluted to a level that does not interfere with the test method and
316 does not exceed the MVD. The PPC is obtained by spiking a test substance with a known EC
317 (e.g., 1 EU/mL) and demonstrating that 50% to 200% of the EC is recovered.

318 **5.6 Negative Product Control (NPC)**

319 The NPC is the test substance diluted to the MVD and then spiked with PFS. It is the
320 negative control for the PPC.

321 **6.0 EXPERIMENTAL DESIGN**

322 **6.1 Incubation with Test Samples and Measurement of IL-6 Release**

323 6.1.1 Preparation of MM6 Cells for the Incubation Plate

324 To perform an ELISA on the following day, obtain 30 to 50 mL of MM6 cell suspension
325 from propagation flasks and centrifuge at 100 x g for 8 min at RT. Remove the supernatant,
326 resuspend the cell pellet in 3-5 mL of RPMI-C (containing 2% FCS for all ELISA
327 procedures), and gently pipet up and down to mix. Transfer aliquots of cells to new culture

328 flasks and add RPMI-C to obtain a concentration of 4×10^5 cells/mL. Place the flasks in a cell
 329 culture incubator and maintain the cells at $37 \pm 1^\circ\text{C} + 5\% \text{CO}_2$. In general, each 96-well assay
 330 plate requires approximately 10 mL of cell suspension at 2.5×10^6 cells/mL.

331 6.1.2 Incubation Plate

332 Test substances are prepared at a level of dilution that did not show interference with the test
 333 system or for which it is known that interference does not occur. Each incubation plate can
 334 accommodate an endotoxin standard curve, a NSC, and 14 test substances (see **Table 6-1**).

335 **Table 6-1 Overview of Incubation Plate Preparation in the MM6/IL-6 Test Method**

Number of Wells	Sample	RPMI-C	EC	Test Sample	MM6 ¹	Mix the samples; incubate overnight at $37 \pm 1^\circ\text{C}$ in a humidified atmosphere with 5% CO_2 .	Mix the samples; immediately transfer to an ELISA plate ⁴ and run ELISA or store plate at -20°C or -80°C .
		μL					
20 ²	EC	100	50	0	100		
4	NSC	150	0	0	100		
56 ³	Test samples (1-14)	100	0	50	100		

336 Abbreviations: EC = Endotoxin control; NSC = Negative saline control; MM6 = Mono Mac 6 cell line

337 ¹ MM6 cell concentration is 2.5×10^6 cells/mL.

338 ² Five EC concentrations (0.125, 0.25, 0.50, 1.0, and 2.0 EU/mL) in quadruplicate.

339 ³ 14 test samples (n=4 each) per plate.

340 ⁴ An IL-6 standard curve is prepared in Columns 11 and 12 on the ELISA plate. Therefore, 80 wells are
 341 available for test samples and controls on the incubation plate.

342

343 6.1.3 Incubation Assay for IL-6 Release

344 **Step 1.** Refer to the suggested incubation plate template presented in **Table 6-2**.

345 **Step 2.** Using a pipetter, transfer 100 μL of RPMI-C into each well.

346 **Step 3.** Transfer 50 μL of test sample into the appropriate wells as indicated in the
 347 template.

348 **Step 4.** Transfer 50 μL of the EC (standard curve) and the NSC controls in
 349 quadruplicate into the appropriate wells according to the template.

350 **Step 5.** Transfer 100 μL of a well-mixed MM6 cell suspension into each well.

351 **Step 6.** Place the covered plate in a tissue culture incubator for 16 to 24 hr at
 352 37±1°C in a humidified atmosphere containing 5% CO₂.

353 **Step 7.** Remove 150 µL of the supernatant from each well, without disrupting the
 354 cells, and transfer to a new 96-well ‘transfer’ plate for the IL-6 ELISA.

355 **Table 6-2 Incubation Plate - Sample and Control Template**

	1	2	3	4	5	6	7	8	9	10	11	12
A	EC ¹ 2.0	EC 2.0	EC 2.0	EC 2.0	TS3	TS3	TS3	TS3	TS11	TS11	Void ³	Void
B	EC 1.0	EC 1.0	EC 1.0	EC 1.0	TS4	TS4	TS4	TS4	TS11	TS11	Void	Void
C	EC 0.50	EC 0.50	EC 0.50	EC 0.50	TS5	TS5	TS5	TS5	TS12	TS12	Void	Void
D	EC 0.25	EC 0.25	EC 0.25	EC 0.25	TS6	TS6	TS6	TS6	TS12	TS12	Void	Void
E	EC 0.125	EC 0.125	EC 0.125	EC 0.125	TS7	TS7	TS7	TS7	TS13	TS13	Void	Void
F	NSC	NSC	NSC	NSC	TS8	TS8	TS8	TS8	TS13	TS13	Void	Void
G	TS1 ²	TS1	TS1	TS1	TS9	TS9	TS9	TS9	TS14	TS14	Void	Void
H	TS2	TS2	TS2	TS2	TS10	TS10	TS10	TS10	TS14	TS14	Void	Void

356 Abbreviations: EC = Endotoxin control; NSC = Negative saline control; TS = Test substance

357 ¹ EC value (e.g., EC 2.0) represents the endotoxin concentration in EU/mL.

358 ² TS number (e.g., TS 1) represents an arbitrary sequence for individual test substances.

359 ³ Columns 11 and 12 are reserved for the IL-6 standard curve on the ELISA plate (see **Table 6-3**).

360
 361

362 6.2 ELISA to Measure IL-6 Release

363 6.2.1 IL-6 Standard Curve

364 An IL-6 standard supplied with the ELISA kit is used. IL-6 standards are typically supplied
365 in lyophilized form and should be reconstituted according to the manufacturer's instructions.
366 The stock solution should be diluted in RPMI-C to the following concentrations: 0, 62.5, 125,
367 250, 500, 1000, 2000, and 4000 pg/mL in volumes of at least 500 μ L. Each well on the
368 ELISA plate will receive 50 μ L of an IL-6 blank or standard.

369 6.2.2 ELISA

370 The manufacturer's instructions provided with the ELISA kit should be followed and a
371 typical experimental design is outlined below. If the user prefers to prepare an in-house
372 ELISA, then appropriate modifications and validation of these changes would be necessary.
373 The ELISA should be carried out at RT and therefore all components must be at RT prior to
374 use. Do *not* thaw frozen specimens by heating them in a water bath. A suggested ELISA
375 plate template is shown in **Table 6-3**, which includes a five-point EC standard curve, a NSC,
376 an eight-point IL-6 standard curve (0 to 4000 pg/mL), and 14 test substances in
377 quadruplicate. The EC standard curve, the NSC, and the test sample supernatants are
378 transferred directly from the incubation plate. The IL-6 standard curve is prepared as
379 described in **Section 6.2.1**. An overview of the ELISA plate preparation is shown in **Table 6-**
380 **4**.

381 **Step 1.** After pipetting up and down very carefully three times (avoid detachment
382 of the adherent MM6 cells) to mix the supernatant, transfer 50 μ L from each well
383 of the Incubation Plate (A1-10; H1-10) to the ELISA plate.

384 **Step 2.** Add 50 μ L of each IL-6 standard (0 to 4000 pg/mL) into the respective
385 wells on the ELISA plate.

386 **Step 3.** Add 200 μ L of the enzyme-labeled detection antibody (neat as supplied, or
387 diluted, if necessary) to each of the wells.

388 **Step 4.** Cover the microtiter plate(s) with adhesive film and incubate for 2 to 3 hr
389 at 20 to 25°C.

390 **Step 5.** Decant and wash each well three times with 300 μ L Buffered Wash
 391 Solution and then rinse three times with deionized water. Place the plates upside
 392 down and tap to remove water.

393 **Step 6.** Add 200 μ L of TMB/Substrate Solution to each well and incubate at RT in
 394 the dark for 15 min. If necessary, decrease the incubation time.

395 **Step 7.** Add 50 μ L of Stop Solution to each well.

396 **Step 8.** Tap the plate gently after the addition of Stop Solution to aid in mixing.

397 **Step 9.** Read the OD₄₅₀ within 15 min of adding the Stop Solution. Measurement
 398 with a reference wavelength of 540 to 590 nm is recommended.

399

400 **Table 6-3 ELISA Plate - Sample and Control Template**

	1	2	3	4	5	6	7	8	9	10	11	12
A	EC ¹ 2.0	EC 2.0	EC 2.0	EC 2.0	TS3	TS3	TS3	TS3	TS11	TS11	IL-6 ³ 0	IL-6 0
B	EC 1.0	EC 1.0	EC 1.0	EC 1.0	TS4	TS4	TS4	TS4	TS11	TS11	IL-6 62.5	IL-6 62.5
C	EC 0.50	EC 0.50	EC 0.50	EC 0.50	TS5	TS5	TS5	TS5	TS12	TS12	IL-6 125	IL-6 125
D	EC 0.25	EC 0.25	EC 0.25	EC 0.25	TS6	TS6	TS6	TS6	TS12	TS12	IL-6 250	IL-6 250
E	EC 0.125	EC 0.125	EC 0.125	EC 0.125	TS7	TS7	TS7	TS7	TS13	TS13	IL-6 500	IL-6 500
F	NSC	NSC	NSC	NSC	TS8	TS8	TS8	TS8	TS13	TS13	IL-6 1000	IL-6 1000
G	TS1 ²	TS1	TS1	TS1	TS9	TS9	TS9	TS9	TS14	TS14	IL-6 2000	IL-6 2000
H	TS2	TS2	TS2	TS2	TS10	TS10	TS10	TS10	TS14	TS14	IL-6 4000	IL-6 4000

401 Abbreviations: EC = Endotoxin control; NSC = Negative saline control; TS = Test substance

402 ¹ EC value (e.g., EC 2.0) represents the endotoxin concentration in EU/mL.

403 ² TS number (e.g., TS1) represents an arbitrary sequence for individual test substances.

404 ³ IL-6 values in columns 11 and 12 are in pg/mL.

405

406

407 **Table 6-4 Overview of ELISA Procedure**

Material transfer from Incubation Plate (μL)	Enzyme-labeled Antibody (μL)	Cover the Incubation Plate and incubate for 2 to 3 hr at 20-25°C.	TMB/Substrate Solution (μL)	Incubate 15 min at RT in dark.	Stop Solution (μL)	Read optical density at 450 nm with a 540-590 nm wavelength reference filter.
50	200		200		50	

408 RT = Room temperature

409

410

411 **7.0 EVALUATION OF TEST RESULTS**412 **7.1 OD Measurements**

413 The OD of each well is obtained by reading the samples in a standard microplate
 414 spectrophotometer (i.e., plate reader) using a visible light wavelength of 450 nm (OD_{450}) with
 415 a 540 to 590 nm reference filter (recommended). OD_{450} values are used to determine assay
 416 acceptability and in the decision criteria for the detection of endotoxin in a test substance (see
 417 **Sections 8.0 and 9.0**).

418 **8.0 CRITERIA FOR AN ACCEPTABLE TEST**

419 Obtain the PPC and the corresponding NPC by interference testing of a test substance in the
 420 presence and absence of a fixed quantity of endotoxin (i.e., 1 EU/mL) in quadruplicate. An
 421 EC (five-point standard curve) and a NSC should be included in each experiment. An IL-6
 422 standard curve should be included in each ELISA as shown in the template presented in
 423 **Table 6-3**. An assay is considered acceptable only if the following minimum criteria are met:

- 424 • The quadratic function of the IL-6 standard curve produces an $r \geq 0.95$ and the
 425 OD_{450} of the blank control is below 0.15.
- 426 • The endotoxin standard curve produces OD_{450} values that ascend in a
 427 sigmoidal concentration response.

428 **9.0 DATA INTERPRETATION/DECISION CRITERIA**

429 **9.1 Decision Criteria for Determination of Pyrogenicity**

430 The *t*-test is used to compare the data of a test sample against the data of the EC (0.5 EU/mL)
431 that is performed in parallel. If this test results in a significant *p*-value (i.e., smaller than 1%),
432 then the sample is considered to be non-pyrogenic, and as pyrogenic if otherwise (Hoffmann
433 et al., 2005), as long as the assay acceptability criteria in **Section 8.0** has been met.

434 **10.0 STUDY REPORT**

435 The test report should include the following information:

436 *Test Substances and Control Substances*

- 437
- 438 • Name of test substance
 - 439 • Purity and composition of the substance or preparation
 - 440 • Physicochemical properties (e.g., physical state, water solubility)
 - 441 • Treatment of the test/control substances prior to testing (e.g., vortexing,
442 sonication, warming, resuspension solvent)

442 *Justification of the In Vitro Test Method and Protocol Used*

443 *Test Method Integrity*

- 444
- 445 • The procedure used to ensure the integrity (i.e., accuracy and reliability) of the
446 test method over time
 - 447 • If the test method employs proprietary components, documentation on the
448 procedure used to ensure their integrity from “lot-to-lot” and over time
 - 449 • The procedures that the user may employ to verify the integrity of the
450 proprietary components

450 *Criteria for an Acceptable Test*

- 451
- 452 • Acceptable concurrent positive control ranges based on historical data
 - Acceptable negative control data

453 *Test Conditions*

- 454 • Cell system used
- 455 • Calibration information for the spectrophotometer used to read the ELISA
- 456 • Details of test procedure used
- 457 • Description of any modifications of the test procedure
- 458 • Reference to historical data of the model
- 459 • Description of evaluation criteria used

460 *Results*

- 461 • Tabulation of data from individual test samples

462 *Description of Other Effects Observed*

463 *Discussion of the Results*

464 *Conclusion*

465 *A Quality Assurance Statement for Good Laboratory Practice (GLP)-Compliant Studies*

- 466 • This statement should indicate all inspections made during the study and the
467 dates any results were reported to the Study Director. This statement should
468 also confirm that the final report reflects the raw data.

469 If GLP-compliant studies are performed, then additional reporting requirements provided in
470 the relevant guidelines (e.g., OECD 1998; EPA 2003a, 2003b; FDA 2003) should be
471 followed.

472

472 **11.0 REFERENCES**

- 473 EPA. 2003a. Good Laboratory Practice Standards. Toxic Substances Control Act. 40 CFR
474 792.
- 475 EPA. 2003b. Good Laboratory Practice Standards. Federal Insecticide, Fungicide, and
476 Rodenticide Act. 40 CFR 160.
- 477 FDA. 2003. Good Laboratory Practices for Nonclinical Laboratory Studies. 21 CFR 58.
- 478 Hoffmann S, Peterbauer A, Schindler S, Fennrich S, Poole S, Mistry Y, Montag-Lessing T,
479 Spreitzer I, Löschner B, van Aalderen M, Bos R, Gommer M, Nibbeling R, Werner-
480 Felmayer G, Loitzl P, Jungi T, Brcic M, Brügger P, Frey E, Bowe G, Casado J, Coecke S, de
481 Lange J, Mogster B, Næss LM, Aaberge IS, Wendel A, Hartung T. 2005. International
482 validation of novel pyrogen tests based on human monocytoid cells. *J Immunol Methods*
483 298:161-173.
- 484 Organization for Economic Cooperation and Development (OECD). 1998. OECD Series on
485 Principle of Good Laboratory Practice and Compliance Monitoring. No. 1. OECD Principles
486 of Good Laboratory Practice (as revised in 1997). Organization for Economic Co-operation
487 and Development (OECD), ENV/MC/CHEM(98)17. Paris: OECD.
- 488 Taktak YS, Selkirk S, Bristow AF, Carpenter A, Ball C, Rafferty B, Poole S. 1991. Assay of
489 pyrogens by interleukin-6 release from monocytic cell lines. *J Pharm Pharmacol* 43:578-582.
- 490

490 **12.0 TERMINOLOGY AND FORMULA**

491 **12.1 Assay Sensitivity (λ)**

492 For an *in vitro* cell-based assay, the variable λ is defined as the lowest statistically significant
493 point on the standard endotoxin concentration-response curve and represents the relative
494 sensitivity of the test method for the detection of endotoxin (i.e., level of detection).

495 **12.2 Endotoxin Control (EC)**

496 The EC is incubated with MM6 cells and serves as the positive control for the experiment.
497 The results should be compared to historical values to insure that it provides a known level of
498 cytokine release relative to the NSC.

499 **12.3 Endotoxin Limit Concentration (ELC)**

500 The ELC is the maximum allowable concentration of endotoxin for a particular product and
501 is expressed in EU per volume (mL) or weight (mg). The ELC is defined by the FDA or
502 specified in the USP². It is calculated as the product of K/M, where:

503 K is the threshold pyrogen dose for parenteral use in rabbits or humans (5.0 EU/kg). At an
504 injection volume of 10 mL/kg, K is equal to 0.5 EU/mL.

505 M is the larger of the rabbit dose or the maximum human dose administered in one hour as
506 defined below and varies with test substance³.

507 **12.4 Maximum Valid Dilution (MVD)**

508 The MVD is the maximum dilution of a test substance that can be tolerated in a test system
509 without exceeding the ELC, if the test substance must be diluted as a result of assay
510 interference. Dilutions beyond the MVD would not be valid for endotoxin detection in the

² ELC values for most marketed pharmaceutical products are provided in the USP or in other pharmacopoeial or regulatory publications (e.g., European Pharmacopoeia, Japanese Pharmacopoeia, Pharmacopoeial guidelines, FDA publications).

³ Values for most marketed pharmaceutical products are provided in the USP or in other pharmacopoeial or regulatory publications (e.g., European Pharmacopoeia, Japanese Pharmacopoeia, pharmacopoeial guidelines, FDA publications).

511 test system. Calculation of the MVD is dependent on whether or not the ELC for a test
512 substance is published. When the ELC is known, the MVD is:

513
$$\text{MVD} = (\text{ELC} \times \text{Product Potency [PP]})/\lambda$$

514 As an example, for "Cyclophosphamide Injection," the ELC is 0.17 EU/mg, PP is 20
515 mg/mL, and the assay sensitivity is 0.1 EU/mL. The calculated MVD would be 34. The test
516 substance can be diluted no more than 1:34 prior to testing.

517 If the ELC is not known, the MVD is:

518
$$\text{MVD} = \text{PP}/\text{Minimum Valid Concentration (MVC)}$$

519 where, $\text{MVC} = (\lambda \times \text{M})/\text{K}$

520 where, M is the maximum human dose

521 As an example, for "Cyclophosphamide Injection," the PP is 20 mg/mL, M is 30
522 mg/kg, and assay sensitivity is 0.1 EU/mL. The calculated MVC is 0.6 mg/mL and the MVD
523 is 33.3. The test substance can be diluted no more than 1:33 in the assay prior to testing.

524 **12.5 Negative Product Control (NPC)**

525 The NPC is a test sample to which PFS is added. The NPC is the baseline for determination
526 of cytokine release relative to the endotoxin-spiked PPC.

527 **12.6 Negative Saline Control (NSC)**

528 The NSC is MM6 cells (in RPMI-C) incubated with PFS (used for dilution of test substance)
529 and is used as the blank.

530 **12.7 Parenteral Threshold Pyrogen Dose (K)**

531 The value K represents the threshold pyrogen dose for parenteral products for rabbits and
532 humans. Based on experimental data, K is fixed at 5.0 EU/kg. For intrathecal products, K is
533 0.2 EU/kg.

534 **12.8 Positive Product Control (PPC)**

535 The PPC is a test substance spiked with the control standard endotoxin (i.e., 0.5 EU/mL or an
536 amount of endotoxin equal to that which produces ½ the maximal increase in OD from the

537 endotoxin standard curve) to insure that the test system is capable of endotoxin detection in
538 the product as diluted in the assay.

539 **12.9 Product Potency (PP)**

540 The concentration for a test substance is the PP typically expressed as $\mu\text{g/mL}$ or mg/mL .

541 **12.10 Rabbit Pyrogen Test Dose or Maximum Human Dose (M)**

542 The variable M represents the rabbit test dose or the maximum human dose in 1 hr. The
543 variable M is expressed in mg/kg and varies with the test substance. For
544 radiopharmaceuticals, M should be adjusted to account for product activity (radioactive
545 decay) at time administration. An average human standard weight of 70 kg is used for the
546 calculation. If a pediatric dose should be used and it is higher than the adult dose, then it
547 should be used in the formula.

548