

## **Appendix B4**

**Draft ICCVAM Recommended Protocol for Future Studies Using the Whole Blood  
(WB)/Interleukin-1 $\beta$  (IL-1 $\beta$ ) Test Method**

***[This Page Intentionally Left Blank]***

1 **Draft ICCVAM Recommended Protocol for Future Studies Using the Whole Blood**  
2 **(WB)/Interleukin-1 $\beta$  (IL-1 $\beta$ ) Test Method**

3

4

**PREFACE**

5

6 This proposed protocol for the detection of pyrogenicity is based on information obtained  
7 from 1) The European Centre for the Validation of Alternative Methods (ECVAM) WB/IL-  
8 1 $\beta$  Background Review Document (BRD) presented in Appendix A of the draft Interagency  
9 Coordinating Committee on the Validation of Alternative Methods (ICCVAM) Pyrogenicity  
10 Test Method BRD, which includes ECVAM Standard Operating Procedures (SOPs) for the  
11 WB/IL-1 $\beta$  test method, and 2) Information provided to the National Toxicology Program  
12 (NTP) Interagency Center for the Validation of Alternative Toxicological Methods  
13 (NICEATM) by Dr. Thomas Hartung, Head of ECVAM. The ECVAM SOPs are based on  
14 the WB/IL-1 $\beta$  methodology first described by Hartung and Wendel (1996). A table of  
15 comparison between the draft ICCVAM recommended protocol and the ECVAM SOPs is  
16 provided. Future studies using the WB/IL-1 $\beta$  test method may include further  
17 characterization of the usefulness or limitations of the assay for regulatory decision-making.  
18 Users should be aware that the proposed test method protocol might be revised based on  
19 additional optimization and/or validation studies. ICCVAM recommends that test method  
20 users routinely consult the ICCVAM/NICEATM website (<http://iccvam.niehs.nih.gov/>) to  
21 ensure that the most current test method protocol is used.

22

22  
23  
24  
25  
26  
27  
28  
29  
30  
31  
32  
33  
34  
35  
36  
37

***[This Page Intentionally Left Blank]***

37 **Table 1 Comparison of Draft ICCVAM Recommended Test Method Protocol with the ECVAM SOP for the Whole**  
 38 **Blood (WB)/Interleukin-1 $\beta$  (IL-1 $\beta$ ) Pyrogen Test Method**

Protocol Component	Draft ICCVAM Protocol	ECVAM SOP	ECVAM Validation SOP
<b>Test Substance</b>	Test neat or at minimal dilution that produces no interference	Same as ICCVAM protocol	Test at MVD
<b>Incubation Plate (number of control or test groups at n=4 each)</b>	NSC (1)	Same as ICCVAM protocol	Same as ICCVAM protocol
	EC (5)	Same as ICCVAM protocol	EC (2)
	TS (14)	TS (25)	EC spikes (3 x 5)
	PPC <sup>1</sup> (0)	Same as ICCVAM protocol	PPC (3)
	NPC <sup>1</sup> (0)	Same as ICCVAM protocol	NPC (3)
	LTAC (0)	LTAC (1)	Same as ICCVAM protocol
<b>ELISA Plate</b>	Includes seven point IL-1 $\beta$ SC and blank in duplicate	Same as ICCVAM protocol	Same as ICCVAM protocol
<b>Decision Criteria for Interference</b>	0.5 x Median OD <sub>450</sub> <sup>2</sup> of 1 EU/mL EC <2x Median OD <sub>450</sub> of 1 EU/mL EC	0.5 x Median OD <sub>450</sub> of 0.5 EU/mL EC <2x Median OD <sub>450</sub> of 0.5 EU/mL EC	Not applicable (tested at MVD)
<b>Assay Acceptability Criteria</b>	Mean OD <sub>450</sub> of 0.5 EU/mL EC $\geq$ 1.6X Mean OD <sub>450</sub> of NSC	Same as ICCVAM protocol	Same as ICCVAM protocol
	Mean OD <sub>450</sub> of PPC is 50% to 200% of 0.5 EU/mL EC	Same as ICCVAM protocol	Not specified
	Mean OD <sub>450</sub> of NSC $\leq$ 0.15	Same as ICCVAM protocol	Same as ICCVAM protocol
	Quadratic function of IL-1 $\beta$ SC $r^2 \geq 0.95$	Same as ICCVAM protocol	Same as ICCVAM protocol
	EC SC produces OD <sub>450</sub> values that ascend in a sigmoidal concentration response	Same as ICCVAM protocol	Same as ICCVAM protocol
<b>Decision Criteria for Pyrogenicity</b>	OD <sub>450</sub> Test sample $\geq$ OD <sub>450</sub> 0.5 EU/mL	Same as ICCVAM protocol	Same as ICCVAM protocol

39 Abbreviations: EC = Endotoxin control; LTAC = Lipoteichoic acid (LTA) control; MVD = Maximum Valid Dilution; NPC = Negative Product Control; NSC =  
 40 Normal saline control; PPC = Positive Product Control; SC = Standard curve; TS = Test substance

41 <sup>1</sup> PPC and NPC are evaluated during the interference test.

42 <sup>2</sup> Median or mean OD<sub>450</sub> values are corrected (i.e., reference filter reading, if applicable, and NSC are subtracted).

43

44

45

46

47

48

49

50

51

52

53

54

55

56

57

58

***[This Page Intentionally Left Blank]***

59

60

61

62

63

## 64 **1.0 PURPOSE AND APPLICABILITY**

65 The purpose of this protocol is to describe the procedures used to evaluate the presence of a  
66 pyrogen (i.e., Gram-negative endotoxin) in a test substance. The presence of Gram-negative  
67 endotoxin is detected by its ability to induce cytokine IL-1 $\beta$  release from monocytoid cells in  
68 human WB. The quantity of IL-1 $\beta$  released is obtained using an enzyme-linked  
69 immunosorbent assay (ELISA) that includes monoclonal or polyclonal antibodies specific for  
70 IL-1 $\beta$ . Release of IL-1 $\beta$  is measured by incubation of WB with test substances or controls  
71 (i.e., positive, negative). The amount of pyrogen present is determined by comparing the  
72 values of endotoxin equivalents produced by WB cells exposed to the test substance to those  
73 exposed to an internationally-harmonized Reference Standard Endotoxin (RSE)<sup>1</sup> or an  
74 equivalent standard expressed in Endotoxin Units (EU)/mL. Based on a rabbit threshold  
75 pyrogen dose of 0.5 EU/mL, which was established in a retrospective evaluation of rabbit  
76 pyrogen test (RPT) data, a test substance is considered pyrogenic if it induces a level of IL-  
77 1 $\beta$  release equal to or greater than 0.5 EU/mL.

78 The focus of this protocol is on the use of the WB/IL-1 $\beta$  test method, specifically for the  
79 detection of Gram-negative endotoxin in parenteral pharmaceuticals. The relevance and  
80 reliability for non-endotoxin pyrogens (e.g., lipoteichoic acid) has not been demonstrated in a  
81 formal validation study.

## 82 **2.0 SAFETY AND OPERATING PRECAUTIONS**

83 All procedures for procurement of eligible blood donors and blood donations should follow  
84 the regulations and procedures set forth by institutional guidelines for utilization of human  
85 substances, which include but are not limited to blood, tissues, and tissue fluids. Standard  
86 laboratory precautions are recommended including the use of laboratory coats, eye  
87 protection, and gloves. If necessary, additional precautions required for specific study  
88 substances or hazardous chemicals will be identified in the Material Safety Data Sheet  
89 (MSDS).

---

<sup>1</sup> RSEs are internationally-harmonized reference standards (e.g., WHO *E. coli* Lipopolysaccharide [LPS] 94/580 [0113:H10:K-]; USP RSE Lot G3E069; FDA Lot EC-6). Equivalent endotoxins include commercially available *E. coli*-derived LPS Control Standard Endotoxin (CSE) or other *E. coli* LPS preparations that have been calibrated with an appropriate RSE.

90 The stop solution used in the ELISA kit is acidic and corrosive and should be handled with  
91 the proper personal protective devices. If this reagent comes into contact with skin or eyes,  
92 wash thoroughly with water. Seek medical attention, if necessary.

93 Tetramethylbenzidine (TMB) solution contains a hydrogen peroxide substrate and 3, 3', 5, 5'-  
94 TMB. This reagent is a strong oxidizing agent and a suspected mutagen. Appropriate  
95 personal protection should be used to prevent bodily contact.

96 Bacterial endotoxin is a toxic agent (i.e., can induce sepsis, shock, vascular damage,  
97 antigenic response) and should be handled with care. Skin cuts should be covered and  
98 appropriate personal protective devices should be worn. In case of contact with endotoxin,  
99 immediately flush eyes or skin with water for at least 15 min. If inhaled, remove the affected  
100 individual from the area and provide oxygen and/or artificial respiration as needed. Skin  
101 absorption, ingestion, or inhalation may produce fever, headache, and hypotension.

## 102 **3.0 MATERIALS, EQUIPMENT, AND SUPPLIES**

### 103 **3.1 Source of Cells**

104 Leukocytes from WB are the source of cells for cytokine production in the WB/IL-1 $\beta$  test  
105 method (Hartung and Wendel, 1996; Schindler et al., 2006). WB is obtained from healthy  
106 human volunteers who have provided their consent according to established institutional  
107 guidelines. Volunteers are expected not to have taken any drugs (e.g., prescription drugs,  
108 recreational drugs, herbal drugs) and to have been free from illness for at least two weeks  
109 prior to donation.

### 110 **3.2 Equipment and Supplies**

111 For all steps in the protocol, excluding the ELISA procedure, the materials that will be in  
112 close contact with samples and/or blood cells (e.g., pipet tips, containers, solutions) should be  
113 sterile and free from detectable pyrogens.

#### 114 3.2.1 Blood Incubation

##### 115 3.2.1.1 *Equipment*

- 116 • Centrifuge



- 117 • Hood; Bio-safety, laminar flow (recommended)
- 118 • Incubator; cell culture (37 $\pm$ 1 $^{\circ}$ C + 5% CO<sub>2</sub>)
- 119 • Pipetter; multichannel (8- or 12-channel)
- 120 • Pipettors; single-channel adjustable (20 and 200  $\mu$ L)
- 121 • Repeating pipetter
- 122 • Vortex mixer
  
- 123 3.2.1.2 *Consumables*
- 124 • Centrifuge tubes; nonpyrogenic, polypropylene (15 and 50 mL)
- 125 • Combitips; repeating pipetter (1.0 and 2.5 mL)
- 126 • Needle set; Sarstedt multifyly, pyrogen-free, 19 mm, 21 gauge for S-Monovette
- 127 • Plates; microtiter, nonpyrogenic, 96-well, polystyrene, tissue culture
- 128 • Reaction tubes; polypropylene (1.5 mL)
- 129 • Reservoirs; fluid
- 130 • Tips; pipetter, sterile, pyrogen-free (20 and 200  $\mu$ L)
- 131 • Tubes; Sarstedt S-Monovette, 7.5 mL, heparinized for blood collection
  
- 132 3.2.2 ELISA
- 133 3.2.2.1 *Equipment*
- 134 • Microplate mixer
- 135 • Microplate reader (450 nm with an optional reference filter in the range of
- 136 600-690 nm)
- 137 • Microplate washer (optional)
- 138 • Multichannel pipetter
  
- 139 3.2.2.2 *Consumables*
- 140 • Container; storage, plastic

- 141 • Deionized water; nonsterile
- 142 • Plates; microtiter, nonpyrogenic, 96-well, polystyrene, tissue culture
- 143 • Pyrogen-free water (PFW)
- 144 • Reservoirs; fluid
- 145 • Tips; pipetter, nonsterile
- 146 • Tubes; polystyrene (12mL)

### 147 3.2.2.3 *ELISA Kit*

148 An ELISA that measures IL-1 $\beta$  release from Cryo WB is used. A variety of IL-1 $\beta$  ELISA  
149 kits are commercially available and the IL-1 $\beta$  ELISA procedure outlined in this protocol is  
150 intended to serve as an example for using an ELISA kit. If the user prefers to prepare an in-  
151 house ELISA, then additional reagents would be required. The IL-1 $\beta$  ELISA should be  
152 calibrated using an IL-1 $\beta$  international reference standard (e.g., WHO 86/680) prior to use.  
153 The IL-1 $\beta$  cytokine assay kits do not provide the RSE or endotoxin equivalent; therefore, this  
154 reagent must be purchased separately. Results obtained using these products are subject to  
155 the assay acceptability and decision criteria described in **Sections 8.0** and **9.0**. IL-1 $\beta$  ELISA  
156 kit components may include the following:

- 157 • ELISA plates coated with anti-human IL-1 $\beta$  capture antibody; monoclonal or  
158 polyclonal
- 159 • Buffered wash solution
- 160 • Dilution buffer
- 161 • Enzyme-labeled detection antibody
- 162 • Human IL-1 $\beta$  reference standard
- 163 • Pyrogen-free saline (PFS)
- 164 • Stop solution
- 165 • TMB/substrate solution

166 **3.3 Chemicals**

- 167 • Endotoxin (e.g., WHO *E. coli* LPS 2nd International Standard 94/580; USP  
168 RSE *E. coli* LPS Lot G3E069; USP RSE *E. coli* Lot G; FDA *E. coli* Lot EC6)

169 **3.4 Solutions**

170 ELISA solutions are listed in **Section 3.2.**

171 **4.0 ASSAY PREPARATION**

172 All test substances, endotoxin, and endotoxin-spiked solutions should be stored at 4°C.

173 **4.1 Endotoxin Standard Curve**

174 An internationally harmonized RSE or equivalent is used to generate the endotoxin standard  
175 curve. The use of any other *E. coli* LPS requires calibration against a RSE using the WB/IL-  
176 1 $\beta$  test method.

177 A standard endotoxin curve consisting of a Normal Saline Control (NSC) and five RSE  
178 concentrations (0.25, 0.50, 1.0, 2.5, and 5.0 EU/mL) are included in the incubation step (refer  
179 to **Table 4-1**) and then transferred to the ELISA plate. To prepare the endotoxin standard  
180 curve, first obtain a 2000 EU/mL stock solution by addition of PFW to the lyophilized  
181 content of the stock vial by following the instructions provided by the manufacturer (e.g., for  
182 a vial containing 10,000 EU, 5 mL of PFW is added). To reconstitute the endotoxin, the stock  
183 vial should be vortexed vigorously for at least 30 min or sonicated in a bath sonicator for 5  
184 min. The stock solution is stable for 14 days when stored at 2 to 8°C. An endotoxin standard  
185 curve is prepared by making serial dilutions of the stock solution in PFS as described in  
186 **Table 4-1.**

187

187

188 **Table 4-1 Preparation of Endotoxin Standard Curve**

Stock Endotoxin EU/mL	$\mu$ L of Stock Endotoxin	$\mu$ L of PFS	Endotoxin Concentration EU/mL
2000 <sup>1,2</sup>	50	1950	50 <sup>3</sup>
50	100	900	5.0
5.0	500	500	2.5
2.5	500	500	1.0
1.0	500	500	0.50
0.50	500	500	0.25
0	0	1000	0

189 Abbreviations: EU = Endotoxin units; PFS = Pyrogen-free saline

190 Each stock tube should be vortexed prior to its use to make the subsequent dilution.

191 <sup>1</sup> A 2000 EU/mL stock solution of endotoxin is prepared according to the manufacturer's instructions.192 <sup>2</sup> The stock solution of USP RSE may be stored in aliquots and kept at -20°C for up to 6 months. Do not store  
193 the endotoxin at -80°C.194 <sup>3</sup> This concentration is not used in the assay.

195

196 **4.2 Test Substances**

197 Liquid test substances should be tested neat or, if interference is detected (see **Section 4.2.1**),  
198 diluted in PFS. Solid test substances should be prepared as solutions in PFS or, if insoluble in  
199 saline, dissolved in dimethylsulfoxide (DMSO) then diluted up to 0.5% (v/v) with PFS,  
200 provided that this concentration of DMSO does not interfere with the assay. The test  
201 substances should be vortexed vigorously for at least 30 min or sonicated in a bath sonicator  
202 for 5 min.

203 **4.2.1 Interference Testing**

204 Interference testing must be carried out on any test sample for which no interference  
205 information is available. The purpose of the interference test is to determine the lowest  
206 dilution (i.e., highest concentration) of a test substance from which an endotoxin spike can be  
207 detected (i.e., based on the decision criteria described in **Section 4.2.1.2**. However, to ensure  
208 a valid test, a test substance should not be diluted beyond its Maximum Valid Dilution  
209 (MVD).

210 For many marketed products, values for the MVD and the Endotoxin Limit Concentration  
211 (ELC) are published in the U.S. Pharmacopeia, the European Pharmacopoeia, and/or Food  
212 and Drug Administration (FDA) guidelines. However if one or both of these values are not

213 available, then calculation of the MVD is dependent on the ELC (see **Section 12.3**). If  
214 unknown, the ELC can be approximated by dividing the maximum hourly dose of the  
215 product by the hourly dose received per patient. For example, if a product is used at an  
216 hourly dose of 100 mg per patient, then the ELC would be 350 EU/100 mg or 3.5 EU/mg.

#### 217 4.2.1.1 *Reference Endotoxin for Spiking Test Substances*

218 The WHO-LPS 94/580 [*E. coli* O113:H10:K-] or equivalent international RSE (e.g., USP  
219 G3E069, FDA EC-6) is recommended for preparation of the endotoxin control (EC). If a  
220 different *E. coli* LPS is used and the potency relative to the RSE is not provided, then each  
221 lot must be calibrated against the RSE in the WB/IL-1 $\beta$  test method. For interference testing,  
222 an endotoxin standard curve (see **Section 4.1**) should be included on each plate.

#### 223 4.2.1.2 *Spiking Test Substances with Endotoxin*

224 For interference testing, non-spiked and endotoxin-spiked test substances are prepared in  
225 microplate wells (n=4 replicates) and an *in vitro* pyrogen test is performed. Either PFS or a  
226 fixed concentration (a concentration selected from the middle of the EC standard curve) of  
227 the RSE (i.e., 1 EU/mL) in PFS is added to the test substance in serial two-fold dilutions. An  
228 illustrative example of endotoxin-spiking solutions is shown in **Table 4-2**. For non-spiked  
229 solutions, 200  $\mu$ L of PFS is added to a well followed by 20  $\mu$ L each of the test substance  
230 (neat or at serial dilution) and 20  $\mu$ L of WB. Endotoxin-spiked solutions are prepared by  
231 adding 180  $\mu$ L of PFS to each well followed by 20  $\mu$ L of the test substance (neat or at serial  
232 dilution) and 20  $\mu$ L of WB. Then, 20  $\mu$ L of a 1 EU/mL solution of endotoxin in PFS is added  
233 and the well contents are mixed (see example presented in **Table 4-2**).

234

234

235 **Table 4-2 Preparation of Endotoxin-spiked and Non-spiked Solutions for**  
 236 **Determination of Test Substance Interference in the Incubation and**  
 237 **ELISA Test Systems**

Sample Addition	Spiked	Non-spiked
	$\mu\text{L}/\text{well}^1$	
Pyrogen-free saline	180	200
Endotoxin spike solution <sup>2</sup>	20	0
Test substance (neat and each serial dilution)	20	20
WB	20	20
Total <sup>3</sup>	240	240

238 Abbreviations: WB = Whole blood

239 <sup>1</sup> n=4 replicates each

240 <sup>2</sup> Endotoxin concentration is 1 EU/mL in PFS.

241 <sup>3</sup> A total volume of 240  $\mu\text{L}$  per well is used for the incubation.

242

243

244 The lowest dilution of the test substance that yields an endotoxin spike recovery of 50% to  
 245 200% in the pyrogen test is determined. The optical density (OD) values of the endotoxin-  
 246 spiked and non-spiked test substances are calibrated against the endotoxin calibration curve.  
 247 The resulting EU value of the non-spiked test substance is subtracted from the corresponding  
 248 EU value of the endotoxin-spiked test substance at each dilution. The % recovery for each  
 249 sample dilution is then determined from the endotoxin spike solution concentration set to  
 250 100%. For example, consider the following interference test results in **Table 4-3**:

251 **Table 4-3 Example of Interference Data Used to Determine Sample Dilution for**  
 252 **Assay**

Sample Dilution	% Recovery of Endotoxin Control
None	25
1:2	49
1:4	90
1:8	110

253

254 Based on these results, the dilution of the test substance used in the *in vitro* pyrogen test  
 255 would be 1:4 (i.e., the lowest dilution between 50% and 200% of the 1 EU/mL EC).

256 4.2.2 Interference with ELISA

257 If the data obtained from the experiment in **Section 4.2.1** suggests the presence of  
258 interference, then a subsequent experiment similar to that described in **Section 4.2.1** would  
259 need to be performed to confirm that the test substance(s) does not directly interfere with the  
260 ELISA. For this experiment, an ELISA would be performed in the absence of PBMCs.

261 **5.0 CONTROLS**

262 **5.1 Negative Control**

263 A negative control (e.g., PFS is added instead of the test sample) is included in each  
264 experiment in order to detect nonspecific changes in the test system, as well as to provide a  
265 baseline for the assay endpoints.

266 **5.2 Solvent Control**

267 Solvent controls are recommended to demonstrate that the solvent is not interfering with the  
268 test system when solvents other than PFS are used to dissolve test substances.

269 **5.3 Positive Control**

270 An internationally standardized EC (e.g., WHO 94/580; USP G3E069; 0.5 EU/mL) is  
271 included in each experiment to verify that an appropriate response is induced.

272 **5.4 Benchmark Controls**

273 Benchmark controls may be used to demonstrate that the test method is functioning properly,  
274 or to evaluate the relative pyrogenic potential of chemicals (e.g., parenteral pharmaceuticals,  
275 medical device eluates) of a specific class or a specific range of responses, or for evaluating  
276 the relative pyrogenic potential of a test substance. Appropriate benchmark controls should  
277 have the following properties:

- 278 • consistent and reliable source(s) for the chemicals (e.g., parenteral  
279 pharmaceuticals, medical device eluates)
- 280 • structural and functional similarities to the class of substance being tested
- 281 • known physical/chemical characteristics

282           • supporting data on known effects in animal models

283           • known potency in the range of response

## 284   **5.5       Positive Product Control (PPC)**

285   The PPC is a test substance diluted to a level that does not interfere with the test method and  
286   does not exceed the MVD. The PPC is obtained by spiking a test substance with a known EC  
287   (e.g., 1 EU/mL) and demonstrating that 50% to 200% of the EC is recovered.

## 288   **5.6       Negative Product Control (NPC)**

289   The NPC is the test substance diluted to the MVD and then spiked with PFS. It is the  
290   negative control for the PPC.

## 291   **6.0       EXPERIMENTAL DESIGN**

### 292   **6.1       Incubation with Test Samples and Measurement of IL-1 $\beta$ Release**

#### 293   6.1.1     Collection of Human Blood

294   WB is obtained from healthy human volunteers who have provided their consent according  
295   to established institutional guidelines. Volunteers are expected not to have taken any drugs  
296   and to have been free from illness for at least two weeks prior to donation. The criteria for  
297   rejection of data from donors that are low responders or that are suspect due to veracity of  
298   health information is addressed in **Section 8.0**.

299   Fresh WB is drawn by venipuncture using a multily needle set and collected in heparinized  
300   tubes (15 IU/mL lithium heparin). All components of the blood collection system (i.e.,  
301   syringes, tubes, connecting lines) must be sterile and pyrogen-free. Blood should be stored at  
302   room temperature (RT) and must be used within 4 hr. Prior to use in the assay, the collection  
303   tubes should be gently inverted 1 to 2 times. **Do Not Vortex**.

#### 304   6.1.2     Incubation Plate

305   Test substances are prepared at a level of dilution that did not show interference with the test  
306   system or for which it is known that interference does not occur. Each incubation plate can  
307   accommodate an endotoxin standard curve, a NSC, and 14 test samples (see **Table 6-1**).

308



308

309 **Table 6-1 Overview of Incubation Plate Preparation in the WB/IL-1 $\beta$  Test Method**

Number of Wells	Sample	PFS	EC	Test Sample	WB	Mix the samples; incubate overnight at 37 $\pm$ 1 $^{\circ}$ C in a humidified atmosphere with 5% CO <sub>2</sub> .	Mix the samples; immediately transfer to an ELISA plate <sup>3</sup> and run ELISA or store plate at -20 $^{\circ}$ C or -80 $^{\circ}$ C.
20 <sup>1</sup>	EC	200	20	0	20		
4	NSC	220	0	0	20		
56 <sup>2</sup>	Test samples (1-14)	200	0	20	20		

310 Abbreviations: EC = Endotoxin control; NSC = Negative saline control; PFS = Pyrogen-free saline; WB =

311 Whole blood

312 <sup>1</sup> Five EC concentrations (0.25, 0.50, 1.0, 2.5, and 5.0 EU/mL) in quadruplicate.313 <sup>2</sup> 14 test samples (n=4) per plate.314 <sup>3</sup> An IL-1 $\beta$  standard curve is prepared in Columns 11 and 12 on the ELISA plate. Therefore, 80 wells are  
315 available for test samples and controls on the incubation plate.

316

317

318 **6.1.3 Incubation Assay for IL-1 $\beta$  Release**

319 Blood samples are prepared in a microtiter plate using a laminar flow hood. All consumables  
320 and solutions must be sterile and pyrogen-free. Each plate should be labeled appropriately  
321 with a permanent marker. An overview of the incubation plate preparation is shown in **Table**  
322 **6-1**. The incubation procedure is outlined below:

323 **Step 1.** Refer to the incubation plate template presented in **Table 6-2**.324 **Step 2.** Using a pipetter, transfer 200  $\mu$ L of PFS into each well.325 **Step 3.** Transfer 20  $\mu$ L of test sample into the appropriate wells as indicated in the  
326 template.327 **Step 4.** Transfer 20  $\mu$ L of the EC (standard curve) and the NSC controls in  
328 quadruplicate into the appropriate wells according to the template.329 **Step 5.** Transfer 20  $\mu$ L of WB into each well and mix by gently swirling the plate.330 **Step 6.** Mix the contents of the wells thoroughly by gently pipetting up and down  
331 five times using a multichannel pipetter, changing the tips between each row to  
332 avoid cross-contamination.

333 **Step 7.** Place the covered plate in a tissue culture incubator for 10 to 24 hr at  
334 37 $\pm$ 1 $^{\circ}$ C in a humidified atmosphere containing 5% CO<sub>2</sub>.

335 **Step 8.** Prior to transferring the test samples onto the ELISA plate, mix the  
336 contents of the wells by pipetting up and down three times using a multichannel  
337 pipetter, changing the tips between each row to avoid cross-contamination.

338 *Note: The aliquots may be tested immediately in the ELISA or stored at -20 $^{\circ}$ C or*  
339 *-80 $^{\circ}$ C for testing at a later time. After transfer to the ELISA plate, freeze the*  
340 *remaining aliquots at -20 $^{\circ}$ C or -80 $^{\circ}$ C for subsequent experiments, if necessary (see*  
341 *Assay Acceptability and Decision Criteria in **Sections 8.0 and 9.0**).*

342 **Table 6-2 Incubation Plate - Sample and Control Template**

	1	2	3	4	5	6	7	8	9	10	11	12
<b>A</b>	EC <sup>1</sup> 5.0	EC 5.0	EC 5.0	EC 5.0	TS3	TS3	TS3	TS3	TS11	TS11	Void <sup>3</sup>	Void
<b>B</b>	EC 2.5	EC 2.5	EC 2.5	EC 2.5	TS4	TS4	TS4	TS4	TS11	TS11	Void	Void
<b>C</b>	EC 1.0	EC 1.0	EC 1.0	EC 1.0	TS5	TS5	TS5	TS5	TS12	TS12	Void	Void
<b>D</b>	EC 0.50	EC 0.50	EC 0.50	EC 0.50	TS6	TS6	TS6	TS6	TS12	TS12	Void	Void
<b>E</b>	EC 0.25	EC 0.25	EC 0.25	EC 0.25	TS7	TS7	TS7	TS7	TS13	TS13	Void	Void
<b>F</b>	NSC	NSC	NSC	NSC	TS8	TS8	TS8	TS8	TS13	TS13	Void	Void
<b>G</b>	TS1 <sup>2</sup>	TS1	TS1	TS1	TS9	TS9	TS9	TS9	TS14	TS14	Void	Void
<b>H</b>	TS2	TS2	TS2	TS2	TS10	TS10	TS10	TS10	TS14	TS14	Void	Void

343 Abbreviations: EC = Endotoxin control; NSC = Negative saline control; TS = Test substance

344 <sup>1</sup> EC value (e.g., EC 5.0) represents the endotoxin concentration in EU/mL.

345 <sup>2</sup> TS number (e.g., TS1) represents an arbitrary sequence for individual test substances.

346 <sup>3</sup> Columns 11 and 12 are reserved for the IL-1 $\beta$  standard curve on the ELISA plate (see **Table 6-3**).

347  
348

## 349 **6.2 ELISA to Measure IL-1 $\beta$ Release**

### 350 **6.2.1 IL-1 $\beta$ Standard Curve**

351 An IL-1 $\beta$  standard, supplied with the ELISA kit, is used. IL-1 $\beta$  standards are typically  
352 supplied in lyophilized form and should be reconstituted according to the manufacturer's  
353 instructions. The stock solution should be diluted in PFS to the following concentrations: 0,

354 62.5, 125, 250, 500, 1000, 2000, and 4000 pg/mL. Each well on the ELISA plate will receive  
355 100  $\mu$ L of an IL-1 $\beta$  blank or standard.

### 356 6.2.2 ELISA

357 The manufacturer's instructions provided with the ELISA kit should be followed and a  
358 typical experimental design is outlined below. If the user prefers to prepare an in-house  
359 ELISA, then appropriate modifications and validation of these changes would be necessary.  
360 The ELISA should be carried out at RT and therefore all components must be at RT prior to  
361 use. Do *not* thaw frozen specimens by heating them in a water bath. A sample ELISA plate  
362 template is shown in **Table 6-3**, which includes a five-point EC standard curve, a NSC, an  
363 eight-point IL-1 $\beta$  standard curve (0 to 4000 pg/mL), and 14 test substances in quadruplicate.  
364 The EC standard curve, the NSC, and the test sample supernatants are transferred directly  
365 from the incubation plate. The IL-1 $\beta$  standard curve is prepared as described in **Section**  
366 **6.2.1**. An overview of the ELISA plate preparation is shown in **Table 6-4**.

367 Immediately prior to the ELISA procedure, dilute or mix any assay components according to  
368 the manufacturer's instructions.

369 **Step 1.** Add 100  $\mu$ L of enzyme-labeled detection antibody to each well.

370 **Step 2.** After pipetting up and down three times to mix the supernatant, transfer  
371 100  $\mu$ L from each well of the Incubation Plate (A1-10; H1-10) to the ELISA plate.

372 **Step 3.** Add 100  $\mu$ L of each IL-1 $\beta$  standard (0 to 4000 pg/mL) into the respective  
373 wells on the ELISA plate.

374 **Step 4.** Cover the microtiter plate(s) with adhesive film and incubate for 90 min on  
375 a microplate mixer at 350-400 rpm at 20 to 25°C.

376 **Step 5.** Decant and wash each well five to six times with 300  $\mu$ L Buffered Wash  
377 Solution per well and then rinse three times with deionized water. Place the plates  
378 upside down and tap to remove the wash solution.

379 **Step 6.** Add 200  $\mu$ L of TMB/Substrate Solution to each well and incubate at RT in  
380 the dark for 10 to 15 min. If necessary, decrease the incubation time.

381 **Step 7.** Add 50  $\mu$ L of Stop Solution to each well.

382 **Step 8.** Tap the plate gently after the addition of Stop Solution to aid in mixing.

383 **Step 9.** Read the OD<sub>450</sub> within 15 min of adding the Stop Solution. Measurement  
384 with a reference wavelength of 600-690 nm is recommended.

385 **Table 6-3 ELISA Plate - Sample and Control Template**

	1	2	3	4	5	6	7	8	9	10	11	12
<b>A</b>	EC <sup>1</sup> 5.0	EC 5.0	EC 5.0	EC 5.0	TS3	TS3	TS3	TS3	TS11	TS11	IL-1 $\beta$ <sup>3</sup> 0	IL-1 $\beta$ 0
<b>B</b>	EC 2.5	EC 2.5	EC 2.5	EC 2.5	TS4	TS4	TS4	TS4	TS11	TS11	IL-1 $\beta$ 62.5	IL-1 $\beta$ 62.5
<b>C</b>	EC 1.0	EC 1.0	EC 1.0	EC 1.0	TS5	TS5	TS5	TS5	TS12	TS12	IL-1 $\beta$ 125	IL-1 $\beta$ 125
<b>D</b>	EC 0.50	EC 0.50	EC 0.50	EC 0.50	TS6	TS6	TS6	TS6	TS12	TS12	IL-1 $\beta$ 250	IL-1 $\beta$ 250
<b>E</b>	EC 0.25	EC 0.25	EC 0.25	EC 0.25	TS7	TS7	TS7	TS7	TS13	TS13	IL-1 $\beta$ 500	IL-1 $\beta$ 500
<b>F</b>	NSC	NSC	NSC	NSC	TS8	TS8	TS8	TS8	TS13	TS13	IL-1 $\beta$ 1000	IL-1 $\beta$ 1000
<b>G</b>	TS1 <sup>2</sup>	TS1	TS1	TS1	TS9	TS9	TS9	TS9	TS14	TS14	IL-1 $\beta$ 2000	IL-1 $\beta$ 2000
<b>H</b>	TS2	TS2	TS2	TS2	TS10	TS10	TS10	TS10	TS14	TS14	IL-1 $\beta$ 4000	IL-1 $\beta$ 4000

386 Abbreviations: EC = Endotoxin control; NSC = Negative saline control; TS = Test substance

387 <sup>1</sup> EC value (e.g., EC 5.0) represents the endotoxin concentration in EU/mL.

388 <sup>2</sup> TS number (e.g., TS1) represents an arbitrary sequence for individual test substances.

389 <sup>3</sup> IL-1 $\beta$  values in columns 11 and 12 are in pg/mL.

390

391

392 **Table 6-4 Overview of ELISA Procedure**

<b>Enzyme-labeled Antibody (<math>\mu</math>L)</b>	<b>Material transfer from Incubation Plate (<math>\mu</math>L)</b>	Incubate 90 min on a plate mixer at 350-400 rpm.	<b>TMB/Substrate Solution (<math>\mu</math>L)</b>	Incubate 15 min at RT in the dark.	<b>Stop Solution (<math>\mu</math>L)</b>	Read optical density at 450 nm with a 600-690 nm wavelength reference filter.
100	100		200		50	

393 Abbreviations: RT = Room temperature

394

395

## 396 7.0 EVALUATION OF TEST RESULTS

### 397 7.1 OD Measurements

398 The OD of each well is obtained by reading the samples in a standard microplate  
399 spectrophotometer (i.e., plate reader) using a visible light wavelength of 450 nm (OD<sub>450</sub>) with  
400 a 600 to 690 nm reference filter (recommended). OD<sub>450</sub> values are used to determine assay  
401 acceptability and in the decision criteria for the detection of endotoxin in a test substance (see  
402 **Sections 8.0** and **9.0**).

### 403 8.0 CRITERIA FOR AN ACCEPTABLE TEST

404 Obtain the PPC and the corresponding NPC by interference testing of a test substance in the  
405 presence and absence of a fixed quantity of endotoxin (i.e., 1 EU/mL) in quadruplicate. An  
406 EC (five-point standard curve) and a NSC should be included in each experiment. An IL-1 $\beta$   
407 standard curve should be included in each ELISA as shown in the template presented in  
408 **Table 6-3**. An assay is considered acceptable only if the following minimum criteria are met:

- 409 • The quadratic function of the IL-1 $\beta$  standard curve produces an  $r^2 \geq 0.95$  and  
410 the OD<sub>450</sub> of the blank control is below 0.15.
- 411 • The endotoxin standard curve produces OD<sub>450</sub> values that ascend in a  
412 sigmoidal concentration response.

### 413 9.0 DATA INTERPRETATION/DECISION CRITERIA

#### 414 9.1 Decision Criteria for Determination of Pyrogenicity

415 The *t*-test is used to compare the data of a test sample against the data of the EC (0.5 EU/mL)  
416 that is performed in parallel. If this test results in a significant *p*-value (i.e., smaller than 1%),  
417 then the sample is considered to be non-pyrogenic, and as pyrogenic if otherwise (Hoffmann  
418 et al., 2005), as long as the assay acceptability criteria in **Section 8.0** has been met.

### 419 10.0 STUDY REPORT

420 The test report should include the following information:

421 *Test Substances and Control Substances*

- 422 • Name of test substance
- 423 • Purity and composition of the substance or preparation
- 424 • Physicochemical properties (e.g., physical state, water solubility)
- 425 • Treatment of the test/control substances prior to testing (e.g., vortexing,
- 426 sonication, warming, resuspension solvent)

427 *Justification of the In Vitro Test Method and Protocol Used*

428 *Test Method Integrity*

- 429 • The procedure used to ensure the integrity (i.e., accuracy and reliability) of the
- 430 test method over time
- 431 • If the test method employs proprietary components, documentation on the
- 432 procedure used to ensure their integrity from “lot-to-lot” and over time
- 433 • The procedures that the user may employ to verify the integrity of the
- 434 proprietary components

435 *Criteria for an Acceptable Test*

- 436 • Acceptable concurrent positive control ranges based on historical data
- 437 • Acceptable negative control data

438 *Test Conditions*

- 439 • Cell system used
- 440 • Calibration information for the spectrophotometer used to read the ELISA
- 441 • Details of test procedure
- 442 • Description of any modifications of the test procedure
- 443 • Reference to historical data of the model
- 444 • Description of evaluation criteria used

445 *Results*

- 446 • Tabulation of data from individual test samples

447 *Description of Other Effects Observed*

448 *Discussion of the Results*

449 *Conclusion*

450 *A Quality Assurance Statement for Good Laboratory Practice (GLP)-Compliant Studies*

451           • This statement should indicate all inspections made during the study and the  
452           dates any results were reported to the Study Director. This statement should  
453           also confirm that the final report reflects the raw data.

454 If GLP-compliant studies are performed, then additional reporting requirements provided in  
455 the relevant guidelines (e.g., OECD 1998; EPA 2003a, 2003b; FDA 2003) should be  
456 followed.

457 **11.0 REFERENCES**

- 458 EPA. 2003a. Good Laboratory Practice Standards. Toxic Substances Control Act. 40 CFR  
459 792.
- 460 EPA. 2003b. Good Laboratory Practice Standards. Federal Insecticide, Fungicide, and  
461 Rodenticide Act. 40 CFR 160.
- 462 FDA. 2003. Good Laboratory Practices for Nonclinical Laboratory Studies. 21 CFR 58.
- 463 Fennrich S, Wendel A, Hartung T. 1999a. New applications of the human whole blood  
464 pyrogen assay (PyroCheck). ALTEX 16:146-149.
- 465 Fennrich S, Fischer M, Hartung T, Lexa P, Montag-Lessing T, Sonntag H-G, Weigand M,  
466 Wendel A. 1999b. Detection of endotoxins and other pyrogens using human whole blood.  
467 Dev Biol Stand. 101:131-139.
- 468 Hartung T, Wendel A. 1996. Detection of pyrogens using human whole blood. *In Vitro*  
469 Toxicol. 9(4):353-359.
- 470 Hoffmann S, Peterbauer A, Schindler S, Fennrich S, Poole S, Mistry Y, Montag-Lessing T,  
471 Spreitzer I, Löschner B, van Aalderen M, Bos R, Gommer M, Nibbeling R, Werner-  
472 Felmayer G, Loitzl P, Jungi T, Brcic M, Brügger P, Frey E, Bowe G, Casado J, Coecke S, de  
473 Lange J, Mogster B, Næss LM, Aaberge IS, Wendel A, Hartung T. 2005. International  
474 validation of novel pyrogen tests based on human monocytoïd cells. *J Immunol Methods*  
475 298:161-173.
- 476 Organization for Economic Cooperation and Development (OECD). 1998. OECD Series on  
477 Principle of Good Laboratory Practice and Compliance Monitoring. No. 1. OECD Principles  
478 of Good Laboratory Practice (as revised in 1997). Organization for Economic Co-operation  
479 and Development (OECD), ENV/MC/CHEM (98)17. Paris: OECD.
- 480 Schindler S, Spreitzer I, Löschner B, Hoffmann S, Hennes K, Halder M, Brügger P, Frey E,  
481 Hartung T, Montag T. 2006. International validation of pyrogen tests based on cryopreserved  
482 human primary blood cells. *J Immunol Methods* 316:42-51.

483

484



484 **12.0 TERMINOLOGY AND FORMULA**

485 **12.1 Assay Sensitivity ( $\lambda$ )**

486 For an *in vitro* cell-based assay, the variable  $\lambda$  is defined as the lowest statistically significant  
487 point on the standard endotoxin concentration-response curve and represents the relative  
488 sensitivity of the test method for the detection of endotoxin (i.e., level of detection).

489 **12.2 Endotoxin Control (EC)**

490 The EC is incubated with WB and serves as the positive control for the experiment. The  
491 results should be compared to historical values to insure that it provides a known level of  
492 cytokine release relative to the NSC.

493 **12.3 Endotoxin Limit Concentration (ELC)**

494 The ELC is the maximum allowable concentration of endotoxin for a particular product and  
495 is expressed in EU per volume (mL) or weight (mg). The ELC is defined by the FDA or  
496 specified in the USP<sup>2</sup>. It is calculated as the product of K/M, where:

497 K is the threshold pyrogen dose for parenteral use in rabbits or humans (5.0 EU/kg). At an  
498 injection volume of 10 mL/kg, K is equal to 0.5 EU/mL.

499 M is the larger of the rabbit dose or the maximum human dose administered in one hour as  
500 defined below and varies with test substance<sup>3</sup>.

501 **12.4 Maximum Valid Dilution (MVD)**

502 The MVD is the maximum dilution of a test substance that can be tolerated in a test system  
503 without exceeding the ELC, if the test substance must be diluted as a result of assay  
504 interference. Dilutions beyond the MVD would not be valid for endotoxin detection in the

---

<sup>2</sup> ELC values for most marketed pharmaceutical products are provided in the USP or in other pharmacopoeial or regulatory publications (e.g., European Pharmacopoeia, Japanese Pharmacopoeia, Pharmacopoeial guidelines, FDA publications).

<sup>3</sup> Values for most marketed pharmaceutical products are provided in the USP or in other pharmacopoeial or regulatory publications (e.g., European Pharmacopoeia, Japanese Pharmacopoeia, Pharmacopoeial guidelines, FDA publications).

505 test system. Calculation of the MVD is dependent on whether or not the ELC for a test  
506 substance is published. When the ELC is known, the MVD is:

507 
$$\text{MVD} = (\text{ELC} \times \text{Product Potency [PP]})/\lambda$$

508 As an example, for "Cyclophosphamide Injection," the ELC is 0.17 EU/mg, PP is 20  
509 mg/mL, and the assay sensitivity is 0.1 EU/mL. The calculated MVD would be 34. The test  
510 substance can be diluted no more than 1:34 prior to testing.

511 If the ELC is not known, the MVD is:

512 
$$\text{MVD} = \text{PP}/\text{Minimum Valid Concentration (MVC)}$$

513 where,  $\text{MVC} = (\lambda \times \text{M})/\text{K}$

514 where, M is the maximum human dose

515 As an example, for "Cylophosphamide Injection," the PP is 20 mg/mL, M is 30  
516 mg/kg, and assay sensitivity is 0.1 EU/mL. The calculated MVC is 0.6 mg/mL and the MVD  
517 is 33.3. The test substance can be diluted no more than 1:33 in the assay prior to testing.

#### 518 **12.5 Negative Product Control (NPC)**

519 The NPC is a test sample to which PFS is added. The NPC is the baseline for determination  
520 of cytokine release relative to the endotoxin-spiked Positive Product Control (PPC).

#### 521 **12.6 Negative Saline Control (NSC)**

522 The NSC is WB and incubated with PFS (used for dilution of test substance) and is used as  
523 the blank.

#### 524 **12.7 Non-intrathecal Threshold Pyrogen Dose (K)**

525 The value K represents the threshold pyrogen dose for parenteral products for rabbits and  
526 humans. Based on experimental data, K is fixed at 5.0 EU/kg. For intrathecal products, K is  
527 0.2 EU/kg.

#### 528 **12.8 Positive Product Control (PPC)**

529 The PPC is a test substance spiked with the control standard endotoxin (i.e., 0.5 EU/mL or an  
530 amount of endotoxin equal to that which produces  $\frac{1}{2}$  the maximal increase in OD from the

531 endotoxin standard curve) to insure that the test system is capable of endotoxin detection in  
532 the product as diluted in the assay.

533 **12.9 Product Potency (PP)**

534 The concentration for a test substance is the PP typically expressed as  $\mu\text{g/mL}$  or  $\text{mg/mL}$ .

535 **12.10 Rabbit Pyrogen Test Dose or Maximum Human Dose (M)**

536 The variable M represents the rabbit test dose or the maximum human dose in 1 hr. The  
537 variable M is expressed in  $\text{mg/kg}$  and varies with the test substance. For  
538 radiopharmaceuticals, M should be adjusted to account for product activity (radioactive  
539 decay) at time administration. An average human standard weight of 70 kg is used for the  
540 calculation. If a pediatric dose should be used and it is higher than the adult dose, then it  
541 should be used in the formula.

542

542

543

544

545

546

547

548

549

550

551

552

553

554

555 ***[This Page Intentionally Left Blank]***

556

557